


AMILens® Individual

Product Data Sheet

Technical Information AMILens Individual (REF: IL01101) / AMILens Individual-Y (REF: ILY1101)

Design:

Acrylic
Custom Aspheric
Sharp Edge



Type

Monofocal, toric posterior chamber IOL
Foldable one piece IOL with 4-loop haptics

Material

Acrylic with UV-absorber; water content 26% at 35°C,
Model *AMILens Individual-Y* with blue light filter

Specifics

Lens to be implanted into the capsular bag of an aphakic eye
Asphericity and dioptric power are customized per eye

Dioptric Range

Customized
Maximum possible power range:
sphere: -20.0 D to +60.0 D
cylinder: +0.5 D to +30.0 D (max. meridional power +60.0 D)

Dioptric Steps

Customized, labelling in steps of 0.1 D

Required biometry data

Corneal tomography, anterior chamber depth, lens thickness and axial length

Manufactured by:

AMIPLANT GmbH, Haidling 1, 91220 Schnaittach, Germany, www.amiplant.com

Phone: +49 9153 924780, E-Mail: mail@amiplant.com