


# AMILens® Individual

## Product Data Sheet

*Technical Information AMILens Individual (REF: IL01101) / AMILens Individual-Y (REF: ILY1101)*

**Design:**

Acrylic  
 Custom Aspheric  
 Sharp Edge



<b>Type</b>	Monofocal, toric posterior chamber IOL Foldable one piece IOL with 4-loop haptics
<b>Material</b>	Acrylic with UV-absorber; water content 26% at 35°C, Model <i>AMILens Individual-Y</i> with blue light filter
<b>Specifics</b>	Lens to be implanted into the capsular bag of an aphakic eye Asphericity and dioptric power are customized per eye
<b>Dioptric Range</b>	Customized Maximum possible power range: sphere: -20.0 D to +60.0 D cylinder: +0.5 D to +30.0 D (max. meridional power +60.0 D)
<b>Dioptric Steps</b>	Labelling in steps of 0.1 D
<b>Required biometry data</b>	Corneal tomography, anterior chamber depth, lens thickness and axial length

**Manufactured by:**

AMIPLANT GmbH, Haidling 1, 91220 Schnaittach, Germany, [www.amiplant.com](http://www.amiplant.com)

Phone: +49 9153 924780, E-Mail: [mail@amiplant.com](mailto:mail@amiplant.com)